English

Valtech Cardio LTD

Cardioband

CARDIOBAND
MITRAL RECONSTRUCTION SYSTEM

INSTRUCTIONS FOR USE

DIRECTORY

1 INTRODUCTION 4

2 INTENDED USE AND INDICATIONS 4

3 CONTRAINDICATIONS 4

4 INTENDED USERS 4

5 LIST OF ACROYNMS 4

6 GENERAL WARNINGS AND CAUTIONS 5

7 POTENTIAL COMPLICATIONS 5

8 PACKAGE CONTENTS 6

9 HOW SUPPLIED 7

10 ADDITIONAL EQUIPMENT NOT PROVIDED WITH SYSTEM 7

11 DEVICE DESCRIPTION 7

12 CARDIOBAND SYSTEM PROCEDURE STEP-BY-STEP INSTRUCTIONS 10

12.1 IMPLANT SIZING 10
12.2 PATIENT PREPARATION 10
12.3 ACCESS TO THE MITRAL VALVE 10
12.4 TSS INSERTION 11
12.5 IMPLANT DELIVERY SYSTEM INSERTION 11
12.6 SYSTEM NAVIGATION ALONG THE POSTERIOR ANNULUS 12
12.6.1 1st, Anchor Placement 13
12.6.2 Consecutive Anchor Deployment 15
12.7 DISENGAGEMENT OF IMPLANT FROM CATHETER 17
12.8 IDS REMOVAL 17
12.9 SAT INSERTION AND ATTACHMENT 18
12.9.1 Implant Adjustment 19
12.10 SAT DISCONNECTION 21
12.11 TSS REMOVAL 21

13 ANTICOAGULANT AND/OR ANTIPLATELET THERAPY 22

14 CLEANING AND STERILIZATION 22

14.1 BOOSTER 22
14.2 CARDIOBAND STAND AND GC-IC FIXATOR 22
14.2.1 Cleaning 22
14.2.2 Steam Sterilization: 23

15 STORAGE 23

16 MAGNETIC RESONANCE (MR) – MR CONDITIONAL 23

17 SPECIFICATIONS 24

17.1 CARDIOBAND SYSTEM DIMENSIONS 24
17.2 CATALOG NUMBERS 25

18 SYMBOL LEGEND 26
19 SYSTEM FIGURES 28

READ ALL INSTRUCTIONS CAREFULLY.

1 INTRODUCTION

The Cardioband Mitral Reconstruction System is a sutureless direct annuloplasty designed for mitral valve repair by a Transvenus transseptal procedure. Cardioband is deployed and fixed along the posterior annulus of the mitral valve, and is then sized to remodel the mitral annulus under echocardiography.

2 INTENDED USE AND INDICATIONS

The Cardioband Mitral Reconstruction System is indicated for the reconstruction and/or remodeling of pathological mitral valves.

3 CONTRAINDICATIONS

1. Active bacterial endocarditis
2. Severe organic lesions with retracted chordae or congenital malformations with lack of valvular tissue
3. Heavily calcified annulus or leaflets
4. Subjects in whom transesophageal echocardiography is contraindicated
5. The device may include traces of Nickel and therefore sensitivity should not be implanted in patients with sensitivity to Nickel
6. Patients who cannot tolerate an anticoagulation/antiplatelet regimen
7. Patients with known severe reaction to contrast agents that cannot be adequately premedicated.

4 INTENDED USERS

Physicians trained to perform invasive endovascular and transseptal procedures who have undergone training for the proper use of the Cardioband system.

5 LIST OF ACROYNMS

AC: Anchor Channel
AD: Anchor Drive
CCW: Counter clockwise
CW: Clockwise
GC: Guide Catheter
GW: Guide Wire
IC: Implant Catheter
IDS: Implant Delivery System
ISE: Implant Securing Element
LA: Left Atrium
RO: Radio-Opaque
SAT: Size Adjusting Tool
TF: Transfemoral
TSS: Transseptal Steerable Sheath

6 GENERAL WARNINGS AND CAUTIONS

- This booklet contains important information on the Cardioband Mitral Reconstruction System. Read all instructions carefully. Failure to follow these instructions, warnings and cautions may lead to damage or patient injury.
- Use of the Cardioband Mitral Reconstruction system is restricted to those physicians trained to perform invasive endovascular and transseptal procedures and those who are trained in the proper use of the Cardioband system
- System preparation is only be performed by qualified personnel in accordance with CBTF System Preparation Manual.
- Medical and/or surgical intervention may be required to correct clinical complications associated with the System. If required, the Cardioband implant can be removed using standard surgical techniques. Perform procedures in a facility with on-site cardiac surgery and immediate access to a cardiac operating room.
- Always inspect the system and its packaging to verify no damage has occurred as a result of shipping and handling and that the sterile barrier has not been compromised. Do NOT use the device if damage is detected
- Do NOT use the system after the Use By date stated on the package label
- The Cardioband Mitral Reconstruction implant, delivery system, and anchor drives are single-use and is intended for single patient use only. Re – use and/or re – sterilization of these system components may result in patient cross contamination, infections, malfunction of the device components and/or other serious injury or death.
- Verify that the Cardioband stand and GC-IC Fixator that are supplied NOT sterile have been sterilized on site before use.
- The System is delicate and should be handled carefully. If excessive force is applied or significant resistance is encountered, the System may be damaged, making it unsafe for use.
- Anatomic considerations may interfere with placement of the Cardioband device. Qualified personnel should ensure anatomic compatibility based on imaging prior to placement of the first anchor.
- Use fluoroscopic and echocardiographic (e.g. trans-esophageal (TEE) and trans-thoracic (TTE)) imaging guidance during implantation and adjustment
- Employ sterile techniques
- To avoid thrombus formation, patients should be heparinized sufficiently to maintain an Activated Clotting Time (ACT) of at least 250 seconds throughout the procedure.
- To minimize the potential of air embolization, ensure flushing and proper de-airing techniques.

7 POTENTIAL COMPLICATIONS

The following are potential risks associated with the Cardioband procedure including potential access complications associated with standard cardiac catheterization for the transfemoral access procedure, and the potential risks of local and/or general anesthesia: Death;

Myocardial infarction; Stroke/ transient ischemic attack; Hemolysis; Heart block; Perforation or damage of vessels or tissue, myocardium or valvular structures; Pericardial effusion/cardiac tamponade; Haematoma; Dyskinesia; Infection including endocarditis and septicaemia; Anaphylactic or cardiogenic shock; Arrhythmia; Emboli (thrombus, air, Cardioband device); Fever; Dizziness, Peripheral ischemia or nerve injury; Hypertension/hypotension; Allergic reaction (contrast, Nickel or chromium, Heparin; Anaesthesia reactions, Atrial septal defect requiring intervention, Multi-system organ failure; radiation; Anemia; Wound dehiscence; Prolonged ventilation, angina; Hypotension; hypertension; peripheral ischemia; vessel spasm; vascular injoury; Hemorrhage; Bleeding ; Deep venous thrombus (DVT), Dyspnea, Edema, pain, Nausea/vomiting, Renal insufficiency or failure, Respiratory failure, Pneumonia; fever; Esophageal irritation or perforation; Conversion to standard valve surgery (emergency or elective); Worsening mitral regurgitation; Worsening heart failure, stenosis; atrioventricular block; failure of the valvular apparatus due to progression of disease; inadequate repair of the valvular structure; detachment/dislodgment (partial or full) of the device; migration or malposition of the device requiring intervention; malfunction of the device; Device degeneration; Device thrombosis requiring intervention; Failure to deliver the device to the intended site; low cardiac output; right heart failure; bleeding complications related to use of anticoagulant therapy.

8 PACKAGE CONTENTS

CARDIOBAND IMPLANTS

One (1) implant per size

CARDIOBAND DELIVERY SYSTEM

One (1) Dilator
One (1) Transseptal Steerable Sheath
One (1) Guide Catheter
One (1) Implant Catheter
One (1) Size Adjustment Tool
One (1) Torque Limiter
One (1) Silicone Mat
One (1) Stand lock screw
One (1) Stand support pin
Two (2) Hemostasis valves
One (1) Stand Cradle
One (1) Wire Feeder
One (1) Peel Away Catheter and Dilator
One (1) Flow Regulator
Five (5) Three-way stopcocks with extension line
Identification Labels
One (1) Instructions for use

CARDIOBAND ANCHOR DRIVES

Four (4) Anchor Drives Kits (16 anchors)

CARDIOBAND SYSTEM ACCESSORIES

One (1) Booster Kit including:
- Booster (1)
- Booster support plate (1)

One (1) Cardioband Stand
One (1) GC-IC Fixator

9 HOW SUPPLIED

Sterile: The Cardioband Implant, Cardioband Delivery System, and Cardioband Anchor drives are provided sterile in sealed pouches. The devices are sterilized with ethylene oxide gas. **Non Sterile**: The Cardioband stand and GC-IC Fixator are provided non-sterile and should be sterilized prior to each use according to "Cleaning and sterilization" procedure. A Booster Kit is provided and should be cleaned after each use.

10 ADDITIONAL EQUIPMENT NOT PROVIDED WITH SYSTEM

Transseptal kit and equipment
Step-up dilators
High pressure extension tubing (4)
1000ml pressure bags (4)
Sterile IV tubing with thumbwheel occluders (4)
1000mL Heparinized Sterile Saline Solution (4)
0.014" High Torque IRONMAN Guidewire, 190cm (1)
Bowl 1000CC (2)
Bowl 500CC (2)
Luer Lock Syringe 20mL (2)

11 DEVICE DESCRIPTION

The Cardioband System consists of the following main parts:

1) Cardioband Implant **(Figure** 1**)**
2) Cardioband Delivery System **(Figure** 2**)** comprised mainly of the following:
 a. Transseptal Steerable Sheath **(Figure** 3**)**
 b. Dilator;
 c. Implant Delivery System, which includes a steerable Guide Catheter **(Figure** 4**)** and an Implant Catheter **(Figure** 5**)**
 d. Size Adjustment Tool **(Figure** 6**)**
3) Cardioband Anchor Drives **(Figure**7**)**

Cardioband Adjustable Implant

Cardioband is a Mitral Reconstruction device designed for mitral valve repair by a Transvenus transseptal procedure. Cardioband is deployed and fixed along the posterior annulus of the mitral valve, and is then adjusted to remodel the mitral annulus.

The Cardioband implant consists of a contraction wire and polyester fabric covering that are attached to an Adjustment mechanism. The polyester covering is marked with radio-opaque (RO) markers at set distances along the implant.

The implant is supplied in several sizes compatible for posterior circumference ranging from 73 to 120 mm.

Cardioband Delivery System

The Cardioband Delivery system is designed to introduce and deploy the Cardioband implant in the appropriate location along the mitral annulus and then adjust the Cardioband implant in order to remodel the annulus.

The Dilator is used for introduction of the TSS to the femoral vein and left atrium. The TSS and GC provide conduits to access the mitral valve and to position the implant on the annulus.

The Implant Catheter enables the user to release portions of the implant fabric along the annulus.

The implant is positioned using anchors. After the implantation, the Implant Delivery system is removed and the Cardioband implant is cinched using the Size Adjustment Tool (SAT). The TSS, Implant Delivery system and SAT are steered and/or actuated by the use of control knobs, handles, latches and levers.

Cardioband Anchor Drives

The anchors are attached to the tip of the Anchor drives. The Anchor Drives are inserted through the Implant Catheter one-by-one, thereby delivering the anchors to the Implant.

Cardioband Stand

The stand **(Figure** 8**)** is positioned over the lower part of the patient's body and is used to support the proximal end of the system tools during the procedure. It allows axial (in/out) movement of the entire system. The Stand is also comprised of two additional components: IC stopper **(Figure** 8**)** that limits the backward motion of the IC and Stand-SAT adaptor **(Figure** 9**)** to enable the placement of the SAT over the stand in the adjustment phase of the procedure. The Stand is provided non-sterile and is sterilized by the facility prior to use.

CardioBand System Handle Controls

Device	Control	Function
Transseptal Steerable Sheath (TSS) **(Figure 10)**	1. Steering knob (CW/CCW) 2. TSS Handle (Right/Left)	1. Tip Steering (CW- Curved – towards Posterior medial commissure, CCW- Straight – towards anterior lateral commissure) 2. TSS Rotation (Right – towards atrium roof, Left – apex)
Delivery Guide Catheter (GC) Handle **(Figure** 11**)**	1. Steering GC knob (CW/CCW) 2. Advancement GC knob	1. Tip steering (into/away from the valve) 2. Advancement and withdrawal maneuvering of the GC tip relatively to the TSS tip
Implant Catheter Handle **(Figure** 12**)**	1. Implant Release knob 2. Implant Catheter Handle 3. Implant Release Safety Latch	1. Releases Implant portions 2. Push/Pull of the Implant Catheter 3. Disconnect implant from Implant Catheter
Anchor Drive **(Figure** 13**)**	1. Anchor drive handle 2. Release lever	1. Anchor Drive rotation 2. Anchor Release
Size Adjustment Tool (SAT) **(Figure** 14**)**	1. Safety latch 2. Adjustment knob 3. Detachment Knob	1. Engagement of the SAT with the implant Adjustment mechanism 2. Contraction/Expansion of the implant 3. Detachment of the SAT and the implant guidewire from the implant

12 CARDIOBAND SYSTEM PROCEDURE STEP-BY-STEP INSTRUCTIONS

12.1 Implant Sizing

The implant size selection is based on preoperative measurements of the mitral annulus size.

Measurements of mitral valve posterior annular circumference commissure to commissure (mm)	Cardioband Implant size
73-80	A
81-88	B
89-96	C
97-104	D
105-112	E
113-120	F

Table 1 - Cardioband Implant Sizes

Warning: Once the appropriate size of implant is selected, ensure that the implant is mounted on the Implant Catheter only by qualified personnel in accordance with CBTF System Preparation Manual

12.2 Patient Preparation

Note: Visually inspect the stand and GC-IC Fixator before each use

1. Prepare the patient per institution's standard practice for transseptal catheterization
2. Prior to draping the patient, place under the patient's leg in the region of the knee, the booster support plate (**Figure 15**).
3. Place the booster (**Figure 15**) on the booster support plate, so that the front edge of the booster will be 90-100 cm from the patient heart.
4. Adjust the height of the booster
5. Prepare the patient for hemodynamic monitoring, per institution's standard practice.

12.3 Access to the Mitral Valve

Warning: Ensure that the Septum is punctured according to anatomic considerations.

1. Following percutaneous puncture of the femoral vein, access the LA using standard trans-venous, transseptal equipment

2. Heparinize the patient. Activated Clotting Time (ACT) should be at a level of at least 250 seconds throughout the whole procedure.

Warning: Failure to administer heparin once transseptal access has been achieved may result in thrombus formation

3. Place a 0.35" extra stiff exchange length guide wire in the left atrium or pulmonary veins.

12.4 TSS Insertion

Warning: Maintain TSS position across the septum throughout the whole procedure.

Warning: Use echocardiography and/or fluoroscopy for guidance and observation.

1. Pre dilate femoral access with a dilator of 18F or wider.
2. Verify the TSS tip is in a straight position
3. Using standard right heart catheterization procedures, insert the TSS with the Dilator through the femoral vein over a high support guide wire that was previously deployed up to the left atrium (LA). Steer the tip of the TSS into the LA under echo guidance. Pull back the tip of the dilator into the TSS, leaving the guidewire in the LA and confirm the tip of the TSS is past the septum.
4. Remove the dilator and Guidewire while aspirating the TSS
5. Connect the connection slot in the handle of the TSS to the Stand

12.5 Implant Delivery System Insertion

Warning: Verify that no air bubbles exist in the funnel lumen.

Warning: Before insertion of the IDS to the TSS, verify that Ironman is completely covered by the funnel

Warning: Use echocardiography and/or fluoroscopy for guidance and observation.

Caution: Avoid overfeeding of the Implant Guidewire throughout the procedure

1. Flush the IDS until saline fills the Implant protection funnel. Verify the funnel hole is facing up
2. Wait a few seconds for air to exit the funnel and carefully insert the IDS into the TSS until the handle of the GC reaches the Stand.
3. De air the GC by allowing back flush through the stopcock.
4. Place the Ironman guidewire inside the left ventricle.
5. Remove the funnel silicon lock
6. Couple the GC to the TSS by inserting the GC pin into the TSS grove and completing 90 degrees CW rotation to lock position.

Warning: Verify that the GC pin is well locked in the end of the TSS grove. This locking mechanism between the GC and TSS provides the system the capability of steering in two perpendicular planes allowing maneuvering of the system tip to the implantation site.

7. Remove the GC- IC fixator by pushing forward the arrows on both sides **(Figure 16).**
8. Open the Blue Hemostatic Valve in order to have IC movement.

Warning: Verify that the Hemostasis valve is open during anchor drive insertion and removal allowing constant backflow to prevent air emboli.

9. Insert the first Anchor Drive with torque limiter through the IDS (Figure**17**).

12.6 System Navigation along the Posterior Annulus

Warning: Do not proceed with the implantation of the first anchor before verifying that imaging adequately supports implantation. In case of inadequate imaging, abort the procedure.

Warning: Do not release the anchor until verifying it is screwed in the tissue by gentle push-pull test under echocardiography and fluoroscopy.

Warning: Do not deploy anchors in the leaflets.

Warning: Avoid damage to adjacent structure during anchor deployment. The proximity of the Circumflex artery proximity should be noted. If proximity is uncertain, consider performing coronary angiography.

Note: The implantation is performed starting from the anterior lateral segment and finishing at the posterior medial segment.

Warning: Do not deploy anchors without the torque limiter.


12.6.1 1st, Anchor Placement

1. Navigate the Implant Catheter tip to the anterior lateral segment using the 3D Echo. The location can be indicated by the Ironman guidewire using fluoroscopy.
2. Gently pull back the Implant guidewire from the GC handle in order to place the adjustment mechanism as proximal as possible from the Implant tip.
3. Position the tip of the Implant Catheter in the first deployment site.
4. Verify tissue contact and the angle between the IC and the annulus plane by using fluoroscopy under perpendicular view and echocardiography under x – plane view.
5. To deploy the first anchor, rotate the torque limiter CW while gently pressing in under imaging guidance (fluoroscopy under perpendicular view and echocardiography under x – plane view).
6. Verify the anchor is firmly fixed in by performing the push – pull test

Note: Anchor unscrewing for re-positioning can be done at any point before anchor drive is removed by rotating the torque limiter in the CCW direction.

7. The below Image sequence should be followed for first anchor deployment:

8. Release the anchor by pulling the release lever (**Figure 13**).
9. Remove ironman guide wire under fluoroscopy guidance.

Warning: Failure to remove ironman might interfere next steps of Cardioband deployment.

10. Introduce the next anchor drive with the torque limiter through the IDS.

12.6.2 Consecutive Anchor Deployment

Warning: Manage the Implant guidewire with minimum slack but no tensioning. .

Warning: At least two anchors must be deployed before the first RO marker.

Warning: The imaging views specified during anchor deployment are critical for procedural success (see point 9 below).

1. Turning the knob of the Implant Catheter (**Figure** 12) releases the implant fabric. Approximately 4 turns (click) of the knob will align the next RO marker with the IC catheter tip.

Warning: Release implant's fabric ALWAYS under fluoroscopy. The scale is merely an estimate and implant release should always be verified under fluoroscopy.

Warning: The RO marker on the implant should not exceed the IC Channel RO marker

2. A scale on the IC handle will indicate the remaining fabric (Implant Release Indication Scale; **Figure** 12). The scale numbering is an estimate of the number of remaining RO markers and therefore equivalent to the remaining anchors to be deployed.

3. Release the implant fabric by turning the Implant Release knob of the Implant Catheter CW until the next RO marker is aligned with the Implant Catheter tip (visible under fluoroscopy using on face view).

Warning: Release implant's fabric ALWAYS under fluoroscopy. The RO marker on the implant should not exceed the IC Channel RO marker.


4. Navigate the system along the annulus to the next anchoring location by using 3D echo. The catheter location should also be verified by fluoroscopy using on face view.
5. Insert the Anchor Drive with torque limiter through the IDS

Note: Verify by imaging that the angle between the catheter and the implant should not exceed 90°.

6. Verify tissue contact and the angle between the IC and the annulus plane by using fluoroscopy under perpendicular view and echocardiography under x – plane view.
7. Under imaging guidance (fluoroscopy under perpendicular view and echocardiography under x – plane view) rotate the torque limiter CW while gently pressing in.
8. Verify the anchor is firmly fixed in by performing the push – pull test.

Note: Anchor unscrewing for re-positioning can be done at any point before anchor drive is removed by rotating the torque limiter in the CCW direction.

9. The below imaging sequence should be followed in each consecutive anchor deployment:

10. Release the anchor by pulling the release lever.
11. Remove the Anchor drive
12. Introduce the next anchor drive with torque limiter.

Note: Anchor positioning is indicated by the RO markers.

13. For deployment of the consecutive anchors, repeat **Steps 3 - 12**.

14. The automatic release of the Red IC Safety latch **(Figure** 12**)** indicates the last RO marker has been reached. This signifies the last anchor deployment location. At this point the Implant Release Knob is locked.

Warning: Do NOT press the Red IC Safety latch before deployment of the last anchor. This may result in implant disengagement.

15. Deploy the last anchor according to steps 4 - 11.

12.7 Disengagement of Implant from Catheter

1. After deployment of the last anchor, to unlock the implant release knob, press down the red IC safety latch.
2. Turn the implant release knob to disengage the implant from the catheter. Disengagement of the implant is indicated by the red IC Safety Latch pop up and the implant release indication scale will reach the bold line mark.

12.8 IDS Removal

Warning: DO NOT remove the IDS before verifying that the TSS is across the septum and inside the left atrium.

1. Un - flex the tip curves of the IDS and TSS catheters

Warning: Verify that the IDS and TSS are un-flexed prior to IDS removal.

2. Verify the Wire Feeder is in the "Inactive" position (**Figure** 18).
3. Insert the Implant Guidewire through the wire feeder base and pull the guidewire out of the wire feeder handle until the wire feeder is positioned over the GC handle.
4. Secure the Wire Feeder by closing the strap on the GC Handle (**Figure** 18)**.**

Warning: Failure to secure the wire feeder in place may result in system malfunction. Ensure alignment between the Wire Feeder and Implant Guidewire port on the GC to avoid kink.

5. Move the Wire Feeder switch to active position (**Figure** 18).
6. Pull up the Plunger (**Figure** 8) of the Stand Cradle to allow the release of the IDS (GC and IC Device) from the Stand Cradle.
7. Remove the IDS carefully, under fluoroscopy guidance, while simultaneously keep Implant Guidewire management with minimum slack and no tensioning using the wire feeder.

Warning: Avoid overfeeding of the Implant Guidewire during IDS removal. Over slack of the Implant Guidewire during removal of the system and SAT insertion may damage the adjustment mechanism.

Warning: Verify that the TSS remains across the septum and inside the left atrium.

8. While retrieving the IDS out of the TSS, once the GC locking segment is exposed, draw back GC funnel from the TSS.

Warning: Failure to draw back the GC funnel while retrieving the IDS will result in bleeding.

9. Remove the IC stopper (**Figure** 8)

12.9 SAT Insertion and Attachment

1. Place the SAT Stand adaptor (located on the stand legs, **Figure** 9) on top of the stand. Loosely fix the SAT Stand adaptor screw.
2. Insert the SAT through the TSS over the Implant guidewire.
3. Once the tip of the blue section of the SAT catheter reaches the TSS hemostatic valve, push the SAT funnel (Figure 19) across the TSS Valve.
4. Continue inserting the SAT until the SAT tip reaches the Implant Adjustment mechanism (visible by fluoroscopy).
5. Verify that the black mark on the implant guide wire is completely across the proximal end of the SAT.
6. Place the SAT over the SAT stand adaptor (Figure 9 and Figure 20)

Warning: Verify using two planes (on face and perpendicular) of fluoroscopy distance of 0.5cm-1cm between the pulley and the implant (distance of at least one spool).

Warning: Verify using two planes (on face and perpendicular) of fluoroscopy that the SAT is not pulling the adjustment mechanism or the implant

7. Once the SAT tip is adjacent to the Implant adjustment mechanism, lock the SAT Stand adaptor.
8. Engage the SAT with the implant by moving the engagement latch on the SAT handle backwards and then forward all the way (Figure 14).

Warning: Verify the SAT latch is advanced all the way forward

12.9.1 Implant Adjustment

Warning: Adjustment should be done under echocardiography and fluoroscopy

1. Release the SAT stand adaptor screw.
2. Rotate the Adjustment roller (**Figure** 14) CW for implant contraction under Echocardiography.
3. The level of adjustment is shown quantitatively by the Adjustment counter window (**Figure** 14) and should not exceed the maximum specified in Table 2.
4. Extend the time between adjustments after counter reached 2.5 for at least 4 minutes.

Warning: Do not cinch the Cardioband more than the maximum level indicated in Table 2

Cardioband Implant size	Maximum Number on the Adjustment counter window
A	3.5
B	4
C	4.5
D	5
E	5.5
F	5.5

Table 2: Cardioband Implant Maximum Adjustment

5. In accordance with the echocardiographic results, if needed, the implant may be expanded (after contraction) by turning the Adjustment Roller CCW. Once adequate results are achieved, if the implant was expanded for sizing, in order to return to the contracted state, ensure that the adjustment roller is turned clockwise 6 clicks prior to SAT removal. Note that this will not influence implant sizing.

Warning: disconnection of the SAT must be done only after turning the adjustment roller clockwise (contracted state)

12.10 SAT Disconnection

1. Before disconnection of the SAT, verify using two planes of fluoroscopy that the SAT is not pulling the adjustment mechanism or the implant.

Warning: Verify using two planes (on face and perpendicular) of fluoroscopy that the SAT is not pulling the adjustment mechanism or the implant

2. Uncover the Disconnection Roller, by sliding the Adjustment roller distally (**Figure 21; arrow a**).
3. Rotate the Disconnection Roller CCW (**Figure** 21**; arrow b**) at least 10 turns, to release the SAT with the implant guide wire from the implant.
4. If the SAT has not disconnected, slide the Adjustment Roller proximally to cover the Disconnection Roller and turn the Adjustment Roller 3-6 clicks CCW. The SAT should now be disconnected.
5. Retract the SAT into the TSS. Draw back the SAT funnel from the TSS proximal end. Remove the SAT.

Warning: Failure to draw back the SAT funnel while retrieving the SAT backwards will result in bleeding

12.11 TSS Removal

1. Steer the tip of the TSS to a straight position and retract it into the right atrium

Warning: Verify that the TSS is un-flexed prior to TSS removal.

2. Unlock the TSS from the Stand by retracting the Securing Element (**Figure** 8).
3. Remove the TSS carefully from the body

13 ANTICOAGULANT AND/OR ANTIPLATELET THERAPY

Short-term anticoagulation therapy be necessary after Cardioband implantation. Anticoagulation and other medical therapy should be used per institutional guidelines.

14 CLEANING AND STERILIZATION

14.1 Booster

The Booster should be cleaned after each use according to institutional guidelines.

14.2 Cardioband Stand and GC-IC Fixator

The Cardioband stand and GC-IC Fixator are provided non-sterile and should be sterilized prior to each use according to the following:

Warning: Product supplied NOT sterile. Steam Sterilization should be performed on site before use

Verify all stand single-use accessories (cradle, pin and screw) are removed and discarded.

Warning: Failure to properly clean and dry the stand and GC-IC Fixator may lead to inadequate sterilization

14.2.1 Cleaning

Pre-Cleaning:
Wash the stand and GC-IC Fixator under a stream of warm water for at least 5 minutes to remove debris. Wipe them with a soft cloth.
Cleaning:
Soak the stand and GC-IC Fixator using an enzymatic cleaner (follow cleaner manufacturer's recommendations). Rinse the device under running water for a minimum of 5 minutes. If debris is still visible – repeat this step.
Drying:
Dry the stand and GC-IC Fixator thoroughly by wiping them with a soft cloth.

14.2.2 Steam Sterilization:

Warning: the sterilization parameters are only valid with sterilization equipment which is properly maintained and calibrated

Warning: The stand and GC-IC Fixator have been validated for use in up to 10 Steam Sterilization Cycles

Use only the steam sterilization method for the stand and GC-IC Fixator. Wrap them in standard sterilization wrap.
To achieve the Sterility Assurance Level (SAL) of 10^{-6}, the minimum recommended temperature and exposure time should be:

- Temperature: 119°-123° at 1.036 Bar
- Cycle Time: 30 minutes

15 STORAGE

Handle with care. The System must be stored between 10 °C-40 °C with relative humidity of 20-80%.

16 MAGNETIC RESONANCE (MR) – MR CONDITIONAL

MR Conditional

Non-clinical testing has demonstrated the Cardioband implant and anchors is MR Conditional. A patient can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3 Tesla or less
- Maximum Spatial gradient magnetic field of 1500-Gauss/cm (extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning (i.e., per pulse sequence).
- Normal operating mode of operation for the MR System

MRI – RELATED HEATING

In non-clinical testing, MRI performed for 15 minutes of MR scanning in a 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) with the Cardioband implant and anchors produced a temperature rise of:

- MR system reported, whole body averaged SAR: 2.9-W/kg

- Calorimetry measured values, whole body averaged SAR: 2.7-W/kg
- Highest temperature change: 2.0°C
- Temperature scaled to whole body averaged SAR of 2-W/kg: 1.4°C

ARTIFACT INFORMATION

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Cardioband implant and anchors. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 20-mm relative to the size and shape of the implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	3,830-mm^2	1,152-mm^2	5,446-mm^2	3,011-mm^2
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

17 SPECIFICATIONS


17.1 Cardioband System Dimensions


Component	Dimension
Implant	
Implant length compatible to mitral valve posterior annular circumference	73-120 mm
Transseptal Steerable Sheath and Dilator	
Dilator working length	900 mm
Transseptal Steerable Sheath proximal shaft diameter	25.5 Fr (8,5mm)
Transseptal Steerable Sheath working length	820 mm
Guide Catheter	
Guide Catheter shaft diameter	20Fr (6,6 mm)
Guide Catheter working length	970 mm
Implant Catheter	
Implant Catheter diameter	13Fr (4,2 mm)
Implant Catheter working length	1280 mm
Size Adjustment Tool	
Size Adjustment Tool working length	1120 mm
Implant Guidewire	
Implant Guidewire diameter	0,4 mm (0,016”)
Implant Guidewire length	2600 mm
Anchor Drive	
Anchor Drive length	1680 mm

17.2 Catalog Numbers

CARDIOBAND IMPLANT	REF no. VSU04003, VSU04004,VSU04005, VSU04006, VSU04007, VSU04008
CARDIOBAND SYSTEM	REF no. VSU04001
CARDIOBAND ANCHOR DRIVES	REF no. VSU04002
CARDIOBAND STAND	REF no.VSU04015
BOOSTER KIT	REF no. VSU04017
GC-IC FIXATOR	REF no.VSU04018


18 SYMBOL LEGEND

Symbol	English
	Sterilized Using ETO
	Caution
	Reference number
	Lot Number
	Manufacturer
	Consult instructions for use
	The System cannot be disposed as unsorted municipal waste. Please contact your local distributor for unit disposal.
	Use by: (date)
	Temperature Limitation
	Do not Reuse
	Do not use if package is damaged
	Do Not Resterilize
	Warning
	Authorized Representative
	Conditions of Storage

	Provided non-sterile
SN	Serial Number
	Relative Humidity
	MR Conditional

 Manufacturer



Valtech Cardio Ltd.
3 Ariel Sharon Ave
Or Yehuda 6037606 Israel
Tel: +972-3-5335959




 Authorized Representative





Meditrial s.r.l.
78 Via Savoia
00198 Roma – Italy
Tel: +39-06-45429780






CE 0344





19 SYSTEM FIGURES



No.	Description	Figure
1	Cardioband Implant	
2	Cardioband Delivery System	1. SAT 2. Silicon mat 3. Wire feeder 4. Torque limiter 5. Guide catheter 6. Implant catheter 7. Peel away catheter & dilator 8. Tray 9. TSS 10. Stand pin & screw 11. Three-way stopcocks 12. Flow regulator 13. Homeostasis valves 14. Stand cradle 15. Dilator 
3	Transseptal Steerable Sheath (TSS)	

No.	Description	Figure
4	Guide Catheter (GC)	
5	Implant Catheter (IC)	
6	Size Adjustment Tool (SAT)	

No.	Description	Figure
7	Anchor Drive	
8	Cardioband Stand	
9	Stand with Stand-SAT Adaptor	
10	TSS Handle	

No.	Description	Figure
11	Guide Catheter (GC) Handle	
12	Implant Catheter Handle	
13	Anchor Drive Handle	
14	SAT Control Handle	
15	Booster Kit	

No.	Description	Figure
16	GC – IC Fixator Removal	
17	Assembled Cardioband System	
18	Wire Feeder Activation (left) Inactive; (right) active	
19	SAT with SAT Funnel	

No.	Description	Figure
20	Assembled SAT	
21	Disconnection Roller Uncovered	

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